



## Donning an Open Flap (Proximal Velcro) CJ Socket

### Overview:

Instructions for donning an Open Flap CJ socket with modified gel liner. It should be noted that the CJ Sail will only be attached with Chicago screws on the lateral side of the socket. This enables the patients to easily engage the Velcro on the socket and liner before tightening the Sail.

### Preparation:

Begin by having the patient in a seated or supine position. Orient the liner so the Proximal pile Velcro modification will line up with the hook Velcro in the deepest part of the Rectus femoris or quadriceps channel in the CJ Socket. Roll the liner onto the patient's limb. **(Figs. 1-3)** Ensure that the liner has been trimmed long enough to prevent the socket and Sail from contacting the patient's skin. **(Fig. 4)** Start by opening the Sail so the hook Velcro is unobstructed. If there is a distal locking mechanism, engage that first and then lay the socket onto the liner. Press down on the socket to join the hook and pile Velcro. **(Fig. 5)** Pull the Sail around the posterior thigh. **(Fig. 6)** If the patient prefers, the medial side of the Sail can be tucked inside the medial edge of the socket. (This is not necessary for the Sail to function properly.) **(Fig. 7)** Tighten the closures until the socket matches the immediate volume of the patient's limb. **(Fig. 8)**

